***Instruction Manual  
Manuale di istruzioni  
Manuel d'instructions  
Manual de instrucciones  
Bedienungsanleitung***

## **JP Recirculating Water Vacuum Pump**

F30620198, F30630198, F30640198

### **General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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#### **This unit must be used for laboratory applications only.**

If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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**This unit has been designed and manufactured in compliance with the following standards:**

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**EN 61010-1:2010**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The use of water jet pumps connected to tap is often impeded by a lack of pressure or by their consume which ranges between 200 and 400 L/h. JP is composed by a plastic water tank in which an impeller driven by an electric motor forces the water into four water jets developing suction. The water is continuously recirculated and its consumption is nil.

For the movement of equipment use the two handles located in front and rear of the tank. To remove the motor from tank use the two lateral handles located on the plastic cover. In the upper part of the instrument a tube holder connects to the internal Venturi ejectors which produce a maximum air flow of 35 l/minute.

The air flow can be regulated from 15 to 35 l/minute by the proper knob. The lateral tube holder can be connected to a plastic tube for sending residual traces of not neutralized gases to a fume hood or to outside environment. The pump is equipped by a handy door for the introduction of water and a level indicator for maximum and minimum. A stopcock located in the rear part of the tank allows the suggested periodical drainage.

It is proposed for the removal and neutralization of corrosive vapours and gases evolved during high temperature mineralizations (Kjeldahl, digestions with nitric, sulfuric or perchloric acid). For this application it is suggested to use with the pump, the SMS fume neutralizer.

As a vacuum source it can be connected to laboratory equipment as rotating evaporators, Buchner, Gooch or membrane filters, gel dryers, etc.

The JP vacuum pump produces at the aspiration tubes a residual pressure of 35 mm Hg with a water temperature of 15 °C.

## 2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- JP recirculating water pump
- Power cable
- Instruction manual

The equipment is to be connected to the right electric supply by the electric cable connected to the rear of the upper part by a panel socket with fuse. Ensure that the electric mains is earthed and that the mains socket is easy to reach.

Fill the tank with tap water through the front door up to the maximum level shown by the external indicator. Turn on the instrument by the lighted main switch.

Air flow and vacuum level are regulated by the proper knob.

**Air flow:** from 15 to 35 l/minute

**Vacuum:** from room pressure (1 bar) to 0.04 bar

Place the instrument at about 15 cm from the walls.

**Attention:** If the unit is used for aspiration of corrosive gases, fumes or liquids or for neutralizing acid gases during Kjeldahl mineralizations adopt the following precautions:

- Use the SMS
- Connect JP to SMS and SMS to the fume cap by using the supplied Viton tube. Viton is a copolymer of vinylidene fluoride and hexafluoropropylene withstanding temperatures up to 200 °C and highly resistant to chemicals (oils, solvents, oxidizing agents, strong alkalies and acids) with good mechanical properties (weathering).

## 3. Operating controls

### 3.1 Absorption

The widely applied Kjeldahl method for protein content evaluation uses a high temperature (up to 420 °C) mineralization step with concentrated sulfuric acid and a catalyst.

During the mineralization corrosive toxic gases such as sulfur dioxide (SO<sub>2</sub>) and sulfur trioxide (SO<sub>3</sub>) are evolved.

JP is able to suck the gases and vapours evolved by up to 20 Kjeldahl digestion tubes. The pump can be connected to SMS fume neutralizer which removes acid gases by three steps (condensing, neutralization and absorption) allowing a correct disposal without emission of toxic substances to the environment (air or water drain).

### 3.2 Time related programs A or B

The efficiency of fume removal equipment SMS is related to the time of residence of fumes deriving from digestion in the various stages. If the digestion test tubes with the samples and reagents are introduced in a heating block when this is hot (normally 420 °C) a large amount of acid fumes is developed in the first minutes. In the following minutes the volume of fumes is reduced.

Thus in the first phase of digestions it is necessary to dispose of a large capacity of gas removal that can be reduced lately for a better yield of the system for gas abatement.

This is the reason of the possibility of selection of an optimal sucking power on the JP related to the different phases of digestion and to the number of digested samples. It is possible to select a work program of JP which changes automatically the air flow during the digestion and in relation to the used heating block.

<b>Program</b>	<b>Digester</b>
A	DK6, DK8, DK6/48, DKL8
B	DK20/26, DK20, DK42/26, DKL12, DKL20, DKL42/26

When the pump is turned on by the general switch a continuous operation is possible with a sucking capacity which is selected by the purposely devised knob. The push-button Set located on the front panel allows the choice of two different time related programs (A or B) with programmed air flows.

The operation of one of the two programs is shown by the lighting of a led.

**Time related program A** - During the first 10 min. of operation the air flow is 50% of the maximum and then lowers to 25%

**Time related program B** - During the first 10 min. of operation the air flow is 100% of the maximum and then lowers to 50%

**NOTE:** times and flows are defined on the basis of sample and reagents amounts described by the methods of VELP Scientifica manuals and for the introduction of digestion test tubes in the heating block when the digestion temperature is almost reached.

## 4. Maintenance

In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

- Replace, periodically the tap water in the tank.
- During operation, if necessary, replace the level of water in the tank up to the maximum level and no more to avoid spillages. A level of water lower than the minimum can compromise pump efficiency.

During long idle periods it is suggested to empty completely the tank and to wash with tap water also the immersed rotating part. In the case of use for aspiration of fumes deriving from Kjeldahl mineralizations with a connection to fume neutralizer Scrubber SMS it is suggested to perform this simple maintenance at the end of each working day.

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

## 5. Technical data

Code	F30620198	F30630198	F30640198
Power supply	230V/50Hz +/- 10%	230V/60Hz +/-10%	115V/60Hz +/-10%
Fuses	2 x T3,15A L 250V	2 x T3,15A L 250V	2 x T5A L 250V
Power	220 W	330 W	220 W
Dimensions (WxHxD)	250x400x370 mm		
Weight	10 kg (empty tank)		
Max. air flow	35 NI/1'		
Flow regulation	from 15 to 35 NI/1'		
Vacuum regulation	from 1 to 0.04 bar		
Residual pressure	35 mm Hg with a water temperature in the tank = 15°C		
Liquid volume	Min. 9 L - Max. 11 L		
Working programs	2		
Enviromental temp. range	from 0 to 40°C		
Max humidity	80%		
Pollution degree EN61010-1	2		
Overvoltage category EN61010-1:2010	II		
Max altitude	Up to 2000 m		
Noise level	75 dBA		

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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and satisfies the essential requirements of the following directives:

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2006/42/EC                      Machinery directive  
 2006/95/EC                      Low voltage directive  
 2004/108/EC                      Electromagnetic compatibility directive

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

## Thank you for having chosen VELP!

Since 1983 VELP has offered to professionals in the sector a range of sophisticated and reliable equipment. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification. Instruments are built according to the International norms IEC 1010-1 and to the rules of the CE mark.

Our product lines:

### Food&Feed Line

DK and DKL Digesters  
JP Recirculating water pump for fumes aspiration  
SMS Scrubber  
UDK Distillation units  
SER Solvent extractors  
FIWE Raw fiber extractors  
Dietary fiber extractors  
Oxitest Oxidation test reactor  
NDA 701 Elemental analyzer  
Consumables

### Environment Line

ECO Thermoreactors  
BOD. Determination systems  
Refrigerated thermostats and incubators  
Flocculators  
Overhead mixer  
Mineralization unit for trace heavy metals determination  
Turbidimeter  
Radiation detector

### Stirring Line

Heating magnetic stirrers  
Vertex digital thermoregulator  
Ultraflat magnetic stirrer  
Magnetic stirrers  
Overhead stirrers  
Heating plates  
Vortex mixers  
Homogenizer

### Other Lab Solution

Open circulating bath  
Recirculating water vacuum pump  
Peristaltic pump

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